

SPONSORED BY JENSEN DENTAL



Artistry with Creation

Capoferri draws on his youthful energy, his passion for beauty, and his eye for symmetry to capture the essence of nature and replicate it at the bench.

This expose of Daniele Capoferri's work demonstrates the demand he puts upon himself to precisely restore the natural smile of each patient he serves.

Capoferri is as meticulous in his ceramic artistry as he is in choosing the ceramic materials that best replicate nature. In each of the cases presented here, he restored smiles using Jensen's Creation porcelain.

PHOTOS AND CASES, DANIELE CAPOFERRI/RESTORATIVE DENTIST, RICHARD KINSELD, DDS

PRECISION AND PASSION

Meet Daniele Capoferri, CDT. At the bench or at play, he combines his quest for perfection with his passion for precision.

Swiss born and educated Daniele Capoferri, CDT is a new generation technician with a passion for excellence. Whether coaching and teaching others, or working with the dental team to create elegant dental solutions, Capoferri settles for nothing less than perfection in his art and in the materials he uses. His artistic ability to replicate the symmetry and precision of nature has earned him membership in the prestigious Oral Design organization and invitations to lecture worldwide to both dentists and technicians on his advanced interdisciplinary techniques.

Capoferri operates Swiss Dental Design laboratory, located in Foster City, Calif. (www.capoferridentaldesign.com; daniele@capoferridentaldesign.com). **lab**



Restoration of a fractured tooth on #8 and old crowns on #7 and #9 demonstrates the close working relationship Capoferri maintains with his clients and the patient during all restorative phases of treatment. Note the harmonious papillary position and the zenith of the facial gingival margin.



A fracture of tooth #8 challenged Capoferri to restore this single central and harmonize it with adjacent teeth and the gingiva.

